



JUST RELEASED!



Liner Test Lung compliance curve is linear and predictable.



Other bag-style lung curves are not predictable or realistic.



Pair with IngMar Medical's Parabolic Resistor Ring for increased range of simulation.

Linear Test Lung™ Essential Ventilator Testing & Training


IngMar Medical's Linear Test Lung provides impressive performance and impressive value. The uniquely designed Linear Test Lung is the only bag-style test lung with the ability to create realistic compliance curves throughout a range of settings, performing more like a real patient than other bag-style lungs. Moreover, the Linear Test Lung is priced significantly below similar bag-style lungs – while providing a higher level of performance.

Features

- Compliances C30, C60, C120 (approx.)
- Resistance Rp20
- Tidal volume 1.0 L (max)
- Autoclavable at 120° C
- Order No. 16 00 100

**For More Information
Call: 412.441.8228**

Respiratory Simulation Specialists



INGMAR MEDICAL

PO Box 10106 Pittsburgh, PA 15232 USA
412.441.8228 Toll free 800.583.9910
www.ingarmed.com