

Parabolic Resistor Ring

A Complete Set of Respiratory Resistors


Increase efficiency of ventilator testing and training with IngMar Medical's Parabolic Resistor Ring. This simple but clever tool provides a complete set of 7 respiratory orifices. Select the resistor you need with a quick rotation of the ring – there is no need to disconnect tubing! Because resistors are integrated into the ring, individual resistors cannot be misplaced.

Never search for a missing resistor again!

- Precision laser-cut orifices – Rp5, 10, 20, 50, 100, 200, 500 in a compact 4" ring
- Meets ISO 10651 and ASTM F1100 ventilator testing standards
- Calibration certificate available for standards testing

**For More Information
Call: 412.441.8228**

Respiratory Simulation Specialists



INGMAR MEDICAL

PO Box 10106 Pittsburgh, PA 15232 USA
412.441.8228 Toll free 800.583.9910
www.ingmarmed.com