



Neonatal Demonstration Lung Model


Resistance Settings




Compliance Settings



Leak Settings





| Part No. | Description |
|-----------|-----------------------------------|
| 1131 | 2.5 mm ET-tube resistor |
| 1132 | 3.5 mm ET-tube resistor |
| 1100 | stopcock manifold |
| 10 10 132 | overdistension bracket |
| 1141 | leak cap (white) |
| 1142 | leak cap (red) |
| 10 10 141 | C1.0 bellows (pair) (clear color) |
| 10 10 200 | Carrying pouch (not pictured) |

WARNING !

Always follow Operating Instructions of your ventilator when checking for safety features and performance before use.

Use of this lung model is not intended to substitute for manufacturer prescribed tests and procedures.