

## *Parabolic Resistor Ring*

### Resistance Settings

-  Rp 5
-  Rp 10
-  Rp 20
-  Rp 50
-  Rp 100
-  Rp 200
-  Rp 500



**Intended use**  
Ventilator testing  
to ISO 10651 and ASTM F1100

**Warning!**  
**NOT FOR PATIENT USE**

## Warning!

Always follow Operating Instructions of your ventilator when checking for safety features and performance before use.

Use of this Parabolic Resistor Ring is not intended to substitute for manufacturer prescribed tests and procedures.

## Ordering Information

<b>Part No.</b>	<b>Description</b>
15 40 100	Parabolic Resistor Ring with 22/15 mm ISO tubing adapter
80 15 482	Calibration certificate (available for standards testing)