

## Parabolic Resistor Ring




# A Complete Set of Respiratory Resistors

Increase efficiency of ventilator testing and training with IngMar Medical's Parabolic Resistor Ring. This simple but clever tool provides a complete set of seven respiratory orifices. Select the resistor you need with a quick rotation of the ring – there is no need to disconnect tubing! Because resistors are integrated into the ring, individual resistors cannot be misplaced.

Never search for a missing resistor again!

- Precision laser-cut orifices – Rp5, 10, 20, 50, 100, 200, 500 in a compact 4" ring
- Meets ISO 10651 and ASTM F1100 ventilator testing standards
- Calibration certificate available for standards testing



**INGMAR MEDICAL**  
Respiratory Simulation Specialists