

Linear Test Lung




Essential Ventilator Testing and Training


IngMar Medical's Linear Test Lung provides impressive performance and impressive value. The uniquely designed Linear Test Lung is the only bag-style test lung with the ability to create realistic compliance curves throughout a range of settings. It performs more like a real patient than other bag-style lungs, maintaining a smooth response to progressive inflation, without any "hiccups" due to bulging-out of the bag. Moreover, the Linear Test Lung is priced significantly below similar bag-style lungs – while providing a higher level of performance.

Specifications

- Compliances 10, 30, and 60 mL/cm H₂O (approx.)
- Resistance cm H₂O/L/S 20
- Tidal volume 1.0 L (max)
- Autoclavable at 120° C



Pair with IngMar Medical's Parabolic Resistor Ring for increased range of simulation.



INGMAR MEDICAL
Respiratory Simulation Specialists

IngMar Medical, Ltd. is ISO 9001:2015 certified.

5940 Baum Boulevard
Pittsburgh, PA 15206 USA
412.441.8228 Toll free 800.583.9910
www.ingarmed.com

© 2018 IngMar Medical, Ltd. | 9-27